



**The TeleRein™ is  
obtainable from:**

**Jennifer Howard  
Kingfishers Int. Ltd.  
e-mail:**

**[j.h.telerein@clear.net.nz](mailto:j.h.telerein@clear.net.nz)**

**1710 Oxford Road  
CUST  
North Canterbury  
NEW ZEALAND**

**Phone: +64 3 3125 655**

**Price on application**

**[www.telerein.co.nz](http://www.telerein.co.nz)**

**One year warranty  
Internationally**


**NZ Patent No: 540119**

**[klphotography.co.nz](http://klphotography.co.nz)**

- **TeleRein™ consists of:**
  - A small Transmitter which is attached by Velcro straps to the bridle.**
  - This is so light the horse does not even feel aware of it.**
- **Connecting Cables are fastened by Velcro to the bridle cheek-pieces. These cables link to Sensors which connect to the rider's own reins and to the bit.**
- **The coach's hand-held receiver.**
- **The kit comes with its own Battery-charger and Operators Manual in a light carry bag**



**KINGFISHERS  
INTERNATIONAL LTD  
Proudly presents  
TeleRein™**



**This scientific new  
Coaching Tool:**

**Informs you about the rider's  
Contact with the horse's mouth**

**Diagnoses whether the  
Contact is light, irregular,  
Heavy, or one-sided**

**Helps to develop the rider's  
Sensitivity and feel and so  
Raises awareness of the horse's  
Welfare**

**Supports the coach's opinion  
with an objective finding**

## TeleRein™

### Provides:

- Accurate and immediate information about the weight and consistency of the Contact between rider and horse on each rein.
- It allows the coach to assess the Contact even when the rider's hand is out of sight.
- It gives objective support to the coach, demonstrating the amount of tension on each rein.
- The message is flashed to the hand-held receiver, 16 x 7 cm (6x3inches).
- Two columns of red lights indicate on a scale of 1 –12 the degree of tension involved on each rein.
- The scale records pressure from 300 grams to 13kg. (10 oz to 27 lbs).

**The New TeleRein is Now Smaller, Lighter, Smarter,**

**More Sensitive, Useful and Practical**



Inventor and Siefried with TeleRein and carry bag



Receiver

## TeleRein™

**Will** function up to a distance of 40 -60 meters from the receiver. (130-190ft approx)

It can be used equally successfully on jumping horses or horses working on the flat.

The new TeleRein has the facility to down load data to a lap-top computer – endless exciting research possibilities.

